



Compounding instructions

Step 1

- Prepare a clean, clear working environment
- Decide on the required volume and concentration of CBD that will be dispensed for your patient's individual needs
- Select the appropriate bottle size
- Use the provided table to draw up the appropriate volume of Cannatrek C100™ Oil (Note this concentrate contains 100mg/ml CBD isolate oil in MCT oil)
- Add to bottle

Step 2

- If wanting to flavour, draw up the appropriate volume of desired flavour
- Add to bottle

Step 3

- Draw up the appropriate volume of MCT oil and add to bottle
- Place SealSafe adapter in place
- Place lid on bottle
- Gently shake bottle

-
- This is a Schedule 4 medicine. Dispense with dispensing label including dosing instructions.
 - CBD medicine should be dispensed in the flatpacked PetCann dispensing box provided, with at least 1 x 1.0ml syringe and consumer information provided.
 - As this is a compounded medication, best before should be no more than six months or in line with the expiry date of any of the raw materials, whichever occurs first.
 - Compounded CBD oil is for ORAL use only. Store below 25°C in a cool dry place, avoiding exposure to heat or sunlight.
-